COLOR SELECTOR



CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

ORDERING INFORMATION

- 1. PICK A TOP COLOR Upholstery for table or laminate for cabinet
- 2. CHOOSE A LAMINATE FRAME COLOR Legs, frame and table backs
- 3. SELECT DOOR & DRAWERLAMINATE COLOR **NOTE:** If one color is selected, it will be used on the entire base
- 4. SELECT DOOR & DRAWER PULLS (See page 22 for complete selection) **NOTE:** If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.





3BG-Burgundy

2SGA-Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days

Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.

3NL-Neutral